



Recombinant Mouse Neuromedin-K receptor (Tacr3)

| | |
|--------------------------|--|
| Product Code | CSB-CF023070MO |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P47937 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Sequence | MASVPTGENWTDGTAGVGSHTGNLSAALGITEWLALQAGNFSSALGLPVTSQ APSQVRAN LTNQFVQPSWRIALWSLAYGLVVAVAVFGNLIVIWILAHKRMRTVTNYFLVNLA FSDAS VAAFNTLVNFIYGVHSEWYFGANYCRFQNFPPITAVFASIYSMTAIAVDRYMAII DPLKP RLSATATKIVIGSIWILAFLLAFPQCLYSKIKVMPGRTL CYVQWPEGPKQHFTYH IIVII LVYCFLLIMGVITYTIVGITLWGGEIPGDTCDKYHEQLKAKRKVVKMMIIVVTF AICWL PYHVYFILTAIYQQLNRWKYIQQVYLASFWLAMSSTMYNPIIYCCLNKRFRAGF KRAFRW CPFIQVSSYDELELKTTRFHPTRQSSLYTVSRMESVTVLYDPSEGDPKSSRK KRAVPRD PSANGC SHREFKSASTTSSFISSPYTSVDEYS |
| Source | in vitro E.coli expression system |
| Target Names | Tacr3 |
| Protein Names | Recommended name: Neuromedin-K receptor Short name= NKR Alternative name(s): NK-3 receptor Short name= NK-3R Neurokinin B receptor Tachykinin receptor 3 |
| Expression Region | 1-452 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |