




Recombinant Human Glucagon receptor (GCGR)

Product Code	CSB-CF009316HU
Abbreviation	Recombinant Human GCGR protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47871
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose
Product Type	Transmembrane Proteins
Immunogen Species	Homo sapiens (Human)
Sequence	AQVMDFLFEKWKLYGDQCHHNLSELLPPPTLVCNRTFDKYSCWPDTPANTTA NISCPWYL PWHHKVQHRFVFKRCGPDGQWVRGPRGQPWRDASQCQMDGEEIEVQKEV AKMYSSFQVMY TVGYSLSLGALLLAILGGLSKLHCTRNAIHANLFASFVLKASSVLVIDGLLRTR YSQK IGDDLVSSTWLS DGAVAGCRVA AVFMQYGIVANYCWLLVEGLYLHNLLGLATL PERSFFS LYLGIGWGAPMLFVVPWAVVKCLFENVQCWTSNDNMGFWWILRFPVFLAILIN FFIFVRI VQLLVAKLRARQMHHTDYKFRLAKSTLTLIPLLGVEHVVFVAFVTDEHAQGTLSR AKLFFD LFLSSFQGLLVAVLYCFLNKEVQSELRRRWHRWRLGKVLWEERNTSNHRASS SPGHGPPS KELQFGRGGGSQDSSAETPLAGGLPRLAESPF
Source	in vitro E.coli expression system
Target Names	GCGR
Protein Names	Recommended name: Glucagon receptor Short name= GL-R
Expression Region	26-477
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.