



Recombinant Pig G-protein coupled receptor 39 (GPR39)

Product Code	CSB-CF009812PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B2ZHY2
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	MASPSRPGNDCSHVIDHSHVPEFEVATWIKITLILLFLVIFVVGILGNSVTIRVTQ VLQK KGYLQKEVTDHMVSLACSDILVFLIGMPVEFYIIWNPLTTPSYTVSCKLHSFLF ETCSY ATLLHVLTLSEFYIAICHPFRYKAMSGPCQVKLLIGFVWVTSTLVALPLLFAMG VEYPL VDVPSHRGLSCNRSRHHSEHPETSNSMSVCTNLSSRWTVFQSSIFGAFIYLV VLVSVAF MCWSMMQALQRSKQGTAAKGQQLQRKSESEESRSARRQTIIFLRLIVVTLA ICWMPNQ IRRMMAAAKPKQDWTKAYFKAYMILLPFSDTFFYLSSVVNPLLYNVSSQQFRS VFAQVLR CRLTLPHANQDKRLRAQAASTMDSARSVHRPLIFLASRSNSARRTDKVFLST SQSESEA KPQSKPQLLNHESPESDSVMKPANPATENGIQEHEV
Source	in vitro E.coli expression system
Target Names	GPR39
Protein Names	Recommended name: G-protein coupled receptor 39
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.