



Recombinant Human Endothelin B receptor-like protein 2 (GPR37L1)

Product Code	CSB-CF009811HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O60883
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>APLHLGRHRAETQEQQSRSKRGTEDEEAKGVQQYVPEEWAEYPRPIHPAGL QPTKPLVAT SPNPGKDGTPDSGQELRGNLTGAPGQRLQIQNPLYPVTESSYSAYAIMLLAL VVFAVGI VGNLSVMCIVWHSYYLKSAWNSILASLALWDFLVLFCLPIVIFNEITKQRLLDG VSCRA VPFMEVSSLGVTTFSLCALGIDRFHVATSTLPKVRPIERCQSILAKLAVIWWGSM TLAVP ELLLWQLAQEPAPTMGTLDSCIMKPSASLPESLYSLVMTYQNARMWWYFGCY FCLPILFT VTCQLVTWRVRGPPGRKSECRASKHEQCESQLNSTVVGLTVVYAFCTLPENV CNIVVAYL STELTRQTLDLLGLINQFSTFFKGAITPVLLLCICRPLGQAFLDCCCCCCEECG GASEA SAANGSDNKLKTEVSSSIYFHKPRESPLLPLGTPC</p>
Source	in vitro E.coli expression system
Target Names	GPR37L1
Protein Names	Recommended name: Endothelin B receptor-like protein 2 Short name= ETBR-LP-2 Alternative name(s): G-protein coupled receptor 37-like 1
Expression Region	26-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.