



Recombinant Bovine Endothelin B receptor-like protein 2 (GPR37L1)

Product Code	CSB-CF618164BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q17QD8
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	VWLQQGGHQPVAQEQPDRSRRGAEREDAKGLQQYVPEGWAEYPRPIRPAA LQPTQPWVAA SPSPDRARATGGSGQEPQGNVTGPPGQRPQVQNPLYPVTERSYGAYAVLLL ALLLFAVGI VGSLAVMCIVWHSYYLKSAWNVLASLALWDFLVFFCLPVVTFHEITKQRLLG AVSCRA VPFVEVSSLGVTTFSLCALGIDRFHVATSTLPKARPIEPCPSILAKLAVIWWGSM TLAAP ELLLWQLVREPSAAGTVDTCIMKPSAHLPELSLVLTQYQARMWWSFGCY FCLPVLFT VTCQLVTWRVRGTPGRKPESRPGPQEPRGARPSSTVAGLAAVHALCALPEN VCNVVAAAYL SAALTRQTLELLGLVTQFSTFFKAALTPLLLLCVSRPLGRAFLDCCCCCGEGC GEGCGG RAAPAGPRAKLHTELSASGYFHKPREAPLLALGTPC
Source	in vitro E.coli expression system
Target Names	GPR37L1
Protein Names	Recommended name: Endothelin B receptor-like protein 2 Short name= ETBR-LP-2 Alternative name(s): G-protein coupled receptor 37-like 1
Expression Region	26-482
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.