



Recombinant Rat Endothelin B receptor-like protein 2 (Gpr3711)

Product Code	CSB-CF870767RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9QYC5
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	AATLSLGGHRAKVQEQSRPRRGTKDEGPKEVQHYPVEEWAIEYKPIHPAGL QPTKPLVA TSPNPKDGDATSESGQELRTNLTGTPSQRLQIQNPLYPVTESSYSAYAVMLLA LVVFAVG IVGNLSVMCIVWHSYYLKSAWNSILASLALWDFLVLFCLPIVIFNEITKQRLLGD VSCR AVPFMEVSSLGVTTFSLCALGIDRFHVATSTLPKVRPIERCQSILAKLAVIIVVGS MMLAV PELLLWQLAQEPTPTMGTVDSCIMKPSADLPESLYSLVMTYQNARMWWYFG CYFCLPILF TVTCQLVTWRVRGPPGRKPECRAGRHEQCESQLNSTVVGLTVVYAFCTLPEN ICNIVVAY LSTELTRQTLDLLGLINQFSTFFKGAITPVLLLCICRPLGQAFLDCCCCCCEEC GGASD SSATVSADSKLKAIEVSSSIYFHKPRESPLLPLGTPC
Source	in vitro E.coli expression system
Target Names	Gpr3711
Protein Names	Recommended name: Endothelin B receptor-like protein 2 Short name= ETBR-LP-2 Alternative name(s): G-protein coupled receptor 37-like 1 G-protein coupled receptor CNS2
Expression Region	25-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.