



# Recombinant Danio rerio Phosphatidylserine synthase 1 (ptdss1)

<b>Product Code</b>	CSB-CF018962DIL
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q803C9
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Sequence</b>	<p>MATFRSQTLKDDVNYRMHFRMINEQQVEDITIQFFYKPHTISLLTVTVLSLM  YFAFTR  DDGDPDSNLRVGLILLVSFFLVISVLAFPNGPFTRPHPAIWRIVFGLSVLYFLFLV  FIIF  LNWDQVKALMFWLDPNLRYAKREADVMEYAVNCHVITWERILSHFDIFAFSHF  WGWGMKA  LLIRSYGLCWTISITWELTELFMHLNFAECWWDQVILDILLCNGGGIWLGM  TVCRFL  EMRTYHWASIKDIHSTTGKIKRAVLQFTPASWTYVRWLDPKSSLQRVMGVYLF  MIIWQLT  ELNTFFLKHIFVFPACHALSWCRILFIGIITAPTQRYYAYLTDQCKRVGTQCW  VFGAI  AFLEALACIKFGQDLFSKTQILYVILWLVLCLAFITFLCLYVMVWYAENYGPRQKS  FSECE DSIYSEAGDSVTECKGEFEIDSTTSCSTRKRRDSGDSRTINGMEK</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	ptdss1
<b>Protein Names</b>	Recommended name: Phosphatidylserine synthase 1 Short name= PSS-1 Short name= PtdSer synthase 1 EC= 2.7.8.29 Alternative name(s): Serine-exchange enzyme I
<b>Expression Region</b>	1-465
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.