



Recombinant Chicken Muscarinic acetylcholine receptor M2 (CHRM2)

Product Code	CSB-CF005382CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P30372
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	MNNSTYINSSSENVIALESPYKTIEVVFIVLVAGSLSLVTIIGNILVMVSIKVNRHL QTV NNYFLFSLACADLIIGIFSMNLYTLTYVIGYWPLGPVVC DLWLALDYVVSNASVM NLLII SFDRYFCVTKPLTYPVKRTTKMAGMMIAAAWVLSFILWAPAILFWQFIVGGRTV PDKDCY IQFFSNPAVTFGTAAAFYLPVIIMTVLYWQISRASKSRIKKGKKEAAQNQDPVS PSLVQ GKIVKPNNNNIPTSSDGLEHNKVQNGKTTGESVMENCVQGEEKDSSNDSTSV SVVPSNTK EDEAAKDASQISASQDHLKVENS KLTCIRIVTKSQKGDCCAPTNTTVEIVGTNG DEKQNS VARKIVKMTKQPAKKKPPPSREKKVTRTILAILLAFIITWTPYNMVLINSFCASCI PGT VWTIGYWLCYINSTINPACYALCNATFKKTFKHELLMCHYKNIGATR
Source	in vitro E.coli expression system
Target Names	CHRM2
Protein Names	Recommended name: Muscarinic acetylcholine receptor M2
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.