



Recombinant Rabbit Alpha-1A adrenergic receptor (ADRA1A)

Product Code	CSB-CF001385RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O02824
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	MVFLSGNASDSSNCTHPPAPVNIKAILLGVILGGLILFGVLGNILVILSVACHRH LHSV THYYIVNLAVADLLLTSTVLPFSAIFEILGYWAFGRVFCNIWAAVDVLCCTASIIS LCVI SIDRYIGVSYPLRYPTIVTQRRGLRALLCVWAFSLVISVGPLFGWRQPAPDDETI CQINE EPGYVLFSAALGSFYVPLTIILAMYCRVYVAKRESRGLKSGLKTDKSDSEQVTL RIHRKN APAGGSGVASAKNKTHFSVRLKFSREKKAAKTLGIVVGCFLCWLPPFLVMP GSFFPD FKPPETVFKIVFWLGYLNSCINPIIYPCSSQEFKAFQNVLKIQLRRKQSSKHA LGYTL HAPSQALEGQHKDMVRIPVGSGETFYKISKTDGVCEWKFFSSMPRGSARITVP KDQSACT TARVRSKSFLQVCCCVGPSTPNPGENHQVPTIKIHTISLSENGEEV
Source	in vitro E.coli expression system
Target Names	ADRA1A
Protein Names	Recommended name: Alpha-1A adrenergic receptor Alternative name(s): Alpha-1A adrenoreceptor Short name= Alpha-1A adrenoceptor Alpha-1C adrenergic receptor
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.