



Recombinant Tyramine receptor tyra-2 (tyra-2)

Product Code	CSB-CF615183CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q19084
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MMSSYVMSPVDETYTLFQILKGSALFLLVLWTIFANSLVFIVLYKNPRLQTVPNL LVGNL AFSDLALGLIVLPLSSVYAIAGEWVFPDALCEVFVSADILCSTASIWNLSIVGLDR YWAI TSPVAYMSKRNRKRTAGIMILSVWISSALISLAPLLGWKQTAQTPNLIYEKNNTVR QCTFL DLPSYTVYSATGSFFIPTLLMFFVYFKIYQAFKHRARQIYRQKLA VSSHVIRKHI ESTI LHEISHVLPTSDEFAKEEEEEEDSESSGQVENGLGNGNDAIIEEDECEDEDSD EKRDDHT SMTTVTATVTGPTEAPYMKREAKISKSVPIEKESAIQKREAKPMRSVMAISYEK VKRHK RKERIYRKSLQRKPKAISAAKERRGVKVLGILGCFTVCWAPFFFTMYVLVQFCK DCSPNA HIEMFITWLGYSNSAMNPIIYTVFNRDYQIALKRLFTSEKKPSSTSRV
Source	in vitro E.coli expression system
Target Names	tyra-2
Protein Names	Recommended name: Tyramine receptor tyra-2
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.