



Recombinant Rat Glucose-dependent insulinotropic receptor (Gpr119)

Product Code	CSB-CF752652RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7TQN8
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>MESSFSFGVILAVLTILIIAVNALVVVAMLLSIYKNDGVGLCFTLNLA VADTLIGVA ISG LVTDQLSSSAQHTQKTLCSLRMAFVTSSAAASVLTVM LIAFD RYLAIKQPLRYF QIMNGL VAGGCIAGLWLISYLIGFLPLGVSIFQQTTYHGPCTFFAVFHPRFVLT LSCAGFF PAVLL FVFFYCDMLKIASVHSQHIRKMEHAGAMVGACRPPRPVNDFKAVRTVSVLIGS FTLSWSP FLITSIVQVACHKCCLYQVLEKYLWLLGVGNSLLNPLIYAYWQREVRQQLCHM ALGLLAD GSTQPQIETLKGKEERKKVGRKTLTYTCDAQTLTYTCDAQTLTYTCDAQTLTYTCDA CDTQTLTY TCDAQTLTYTCDAQTLTYTCDAQTLTYTCDAQTLTYTCDAQTLTYTCDAQ TLTYTCDA QTLTYTCDAQTLTYTCDAQTLTYTSSLVTGQTEQTPLKRANMSDPLRTRCG</p>
Source	in vitro E.coli expression system
Target Names	Gpr119
Protein Names	Recommended name: Glucose-dependent insulinotropic receptor Alternative name(s): G-protein coupled receptor 119
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.