



Recombinant *Didelphis marsupialis virginiana* Alpha-2C adrenergic receptor (ADRA2C)

Product Code	CSB-CF001390DKL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P35405
Product Type	Transmembrane Protein
Immunogen Species	<i>Didelphis virginiana</i> (North American opossum) (<i>Didelphis marsupialis virginiana</i>)
Sequence	MDLQLTTNSTDSGDRGGSSNESLQRQPPSQYSPAEVAGLAAVVSFLIVFTIVG NVLVIA VLTSRALKAPQNLQVSLASADILVATLVMPFSLANELMNYWYFGKVWCVIYLA LDVLF TSSIVHLCAISLDRYWSVTQAVEYNLKRTPRRIKGIIVTVWLISAVISFPPLISLYR DPE DDLYPQCELNDETWYILSSCIGSFFAPCIIMVLVYVRIYRVAKLRTRLSEKRTV PEGSS QTENGLSRPPVGAGPSTAAAAAASLRLQAGENGHYHLHHHHHHHHHHHHHHHH HHQLRKS LEDIELESSTSENRRRRRSREEAAARKGSRGFSFSSTKGGQSAGAGSRL SRASNRL EFFSTHRRRKRSSLCRRKVTQAREKRFTFVLAVVMGVFVVCWFPPFFFTYSLY GICREACQ VPETLFKFFFWIGYCNSSLNPVIYTIFNQDFRRSFKHILFKKKKTSLQ
Source	in vitro E.coli expression system
Target Names	ADRA2C
Protein Names	Recommended name: Alpha-2C adrenergic receptor Alternative name(s): Alpha-2C adrenoreceptor Short name= Alpha-2C adrenoceptor Short name= Alpha-2CAR
Expression Region	1-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.