



Recombinant Human Probable G-protein coupled receptor 152 (GPR152)

Product Code	CSB-CF848412HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8TDT2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDTTMEADLGATGHRPRTTELDDEDSYPPGGWDTVFLVALLLLGLPANGLMALMA WLAGSQARH GAGTRLALLLSLALSDFLFLAAAFQILEIRHGGHWPLGTAACRFYYFLWGVSYSSGLF LLAALSIDRCLLALCPHWYPGHRPVRLPLWVCAGVWVLATLFSVPWLVPFPA AVWWYDLV ICLDFWDSEELSLRMLEVLGGFLPFLLLLCHVLTQATACRTCHRQQQPAACRGFARVAR TILSAYVVLRLPYQLAQLLYLAFLWDVYSGYLLWEALVYSDYLILLNSCLSPFLCLMASA DLRTLRLSVLSSFAAALCEERPGSFTPTPTQQLDSEGPTLPEPMAEAQSQMDPVAQPQV NPTLQPRSDPTAQPQLNPTAQPQSDPTAQPQLNLMAQPQSDSVAQPQADTNVQTPAPAAS SVPSPCDEASPTPSSHPTPGALEDPATPPASEGESPSSTPPEAAPGAGPT
Source	in vitro E.coli expression system
Target Names	GPR152
Protein Names	Recommended name: Probable G-protein coupled receptor 152 Alternative name(s): G-protein coupled receptor PGR5
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.