



Recombinant *Drosophila melanogaster* Probable G-protein coupled receptor Mth-like 8 (mthl8)

Product Code	CSB-CF893357DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9W0V7
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	FHEETHYPCAFIDTANITGSYGLDGPVHNWTVIPRHFVAVYDFVIENGIRIPAS RHLRA CVCKTKPCVRICCLRGEIYDLEKRQCLVPVAGVSSLPSSHMEVELGNGSLRL VKLQPRF SIHVETPCEHMKAVTKGSEYVHWTLHENG TISHRGHIFSKHYCFTPLLHGNST WEWQPLA CAPEKLYFVLGVREWTYAICLLIAILSMFIVLMVYLMCSEMRNSFYGVAIKAYAIC MILG YALLAYLTLHNPANLSNAACRILPSLALMNLVLSFYILSFI AFKLYLSFYGVVFTKL MFW LIFTPIVLVAVGWSFFVGFSSYYGSRLIFGGDTCWFDPRNWSVMIYFYAPV FVAC AISGFF YVLSQIYIRDQPDIEKSFESIEKNRFKSFWKYFGYTAVVWVVCICSFAFNYY WENRSH LNYAVSFCMAFHGFAALYALIGKNQQIQNFLRRIDNGEDTCENSVPLSSFG
Source	in vitro E.coli expression system
Target Names	mthl8
Protein Names	Recommended name: Probable G-protein coupled receptor Mth-like 8 Alternative name(s): Protein methuselah-like 8
Expression Region	22-492
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.