



Recombinant Mouse Relaxin-3 receptor 1 (Rxfp3)

Product Code	CSB-CF804323MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8BGE9
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<pre>MQVASATPAATVRKAAAGDELSEFFALTPDLLEVANASGNASLQLQDLWWEL GLELPDGA APGHPPGGGGAESTDTEARVRILISAVYWVVCALGLAGNLLVLYLMKSKQGW RKSSINLF VTNLALTDQFVLTLPFWAVENALDFKWPFKAMCKIVSMVTSMNMYASVFFL TAMSVAR YHSVASALKSHRTRGRGRGDCCGQSLRESCCFSAKVLGGLIWASAALASLPN AIFSTTIR VLGEELCLMHFPDKLLGWDRQFWLGLYHLQKVLLGFLPLSIIISLCYLLLVRFIS DRRVV GTTDAVGAAAAPGGGLSTASARRRSKVTKSVTIVVLSFFLCWLPNQALTTWSI LIKFNAV PFSQEYFQCQVYAFPVSVCLAHSNSCLNPILYCLVRREFRKALKNLLWRIASPS LTNMRP FTATTKPEPEDHGLQALAPLNAAAEPDLIYPPGVVVYSGGRYDLLPSSSAY</pre>
Source	in vitro E.coli expression system
Target Names	Rxfp3
Protein Names	Recommended name: Relaxin-3 receptor 1 Short name= RLN3 receptor 1 Alternative name(s): G protein-coupled receptor SALPR homolog Relaxin family peptide receptor 3
Expression Region	1-472
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.