



Recombinant Bovine Pituitary adenylate cyclase-activating polypeptide type I receptor (ADCYAP1R1)

Product Code	CSB-CF642772BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q29627
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MHSDCIFKKEQAMCLEKIQRVNDLMGLNDSSPGCPGMWDNITCWKPAHVGE MVLVSCPEL FRIFNPDQVWETETIGEFGFADSKSLDLSDMRVVS RNCTEDGWSEPFPHYFD ACGFEEYE SETGDQDYYYLSVKALYTVGYSTSLVTLTTAMVILCRFRKLHCTRNFHIMNLFV SFMLRA ISVFIKDWILYAEQDSNHCFVSTVECKAVMVFFHYCVVS NYFWLFIEGLYLFTLL VETFF PERRYFYWYIIIGWGTP TVCVSVWAMLRLYFDDTGCWDMNDNTALWWWIKGP VVG SIMVN FVLFIGIIVILVQKLQSPDMGGNESSIYFSCVQKCYCKPQRAQQHSCKMSELSTI TLRLA RSTLLIPLFGIHYTVFAFSPENVSKRERLVFELGLGSFQGFVAVLYCFLNGEV QAEIK RKWRSWKVNRYFTMDFKHRHPSLASSGVNGGTQLSILSKSSSQIRMSG LPAD NLAT
Source	in vitro E.coli expression system
Target Names	ADCYAP1R1
Protein Names	Recommended name: Pituitary adenylate cyclase-activating polypeptide type I receptor Short name= PACAP type I receptor Short name= PACAP-R-1 Short name= PACAP-R1
Expression Region	38-513
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.