



Recombinant *Macaca fascicularis* C3a anaphylatoxin chemotactic receptor (C3AR1)

Product Code	CSB-CF764894MOV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6TAC8
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Sequence	MAPFSAETNSTDLLSQPWNESPVILSMVILSLTFLGLPGNGLVLWVAGLKMQR RTVNTVW FLHLTLADLLCCLSLPFSLAHLALQGQWPYGRFLCKLIPSIIVLNMFASVFLITAI TLDR CLVVFKPIWCQNHRNVGTACTICGCIWVVAFMCIPVFVYREIFTTDNHSKCGY KFDLSG SLDYLDFYGNPLENRSLENIFQLPGEMNDRLDPSSFQANDHPWTVTTAFHTQT FQRPSAD SLPKDSASHSPYSNVFKPAGVASPKIPSGFPIEGHKTSPLDDSDAFLSTHLKLF PSAASN SLYEYELPQDFQDYYSLGQFTYDNQVSTPLVAITITRLVVGFLLPSVIMIACYSFI VLRM QRGRFAKSQGKTLRVAVVVVTVFLVCWAPYHISGVLSLFTDPETPLGKTLMSW DHVFTAL ASANSCFNPFLYALLGKDFRKKARQSIQSILEAAFSEELTHSTHCSSNNVFSER NSISTT V
Source	in vitro E.coli expression system
Target Names	C3AR1
Protein Names	Recommended name: C3a anaphylatoxin chemotactic receptor Short name= C3AR Short name= C3a-R
Expression Region	1-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.