



Recombinant Drosophila melanogaster Odorant receptor coreceptor (Orco)

Product Code	CSB-CF894717DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VNB5
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	<p>MTTSMQPSKYTGLVADLMPNIRAMKYSGLFMHNFTGGSAFMKKVYSSVHLVF LLMQFTFI LVNMALNAEEVNELSGNTITTLFFTHCITKFIYLAQNQKNFYRTLNIWNQVNTHP LFAES DARYHSIALAKMRKLFVLMLTTVASATAWTTITFFGDSVKMVDHETNSSIPV EIPRLP IKSFYPWNASHGMFYMISFAFQIYYVLFSMIHSNLCDVMFCSWLIFACEQLQHL KGIMKP LMELSASLDTYRPNAAALFRSLSANSKSELIHNEEKDPGTDMDMSGIYSSKAD WGAQFRA PSTLQSFGGNGGGGNGLVNGANPNGLTKKQEMMVRSIAIKYWVERHKHVRL VAAIGDTYG AALLHMLTSTIKLTLLAYQATKINGVNVYAFTVVGYLGYALAQVFHFCIFGNRLI EESS SVMEAAYSCHWYDYGSEEAKTFVQIVCQQCQKAMSISGAKFFFTVSLDLFASVL GAVVTYFM VLVQLK</p>
Source	in vitro E.coli expression system
Target Names	Orco
Protein Names	Recommended name: Odorant receptor coreceptor Alternative name(s): Odorant receptor 83b
Expression Region	1-486
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.