



Recombinant Rat Histamine H1 receptor (Hrh1)

Product Code	CSB-CF010737RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P31390
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>MSFANTSSTFEDKMCEGNRTAMASPQLLPLVVVLSSISLVTVGLNLLVLYAVHS ERKLHT VGNLYIVSLSVADLIVGAVVMPMNILYLIMTKWLSLGRPLCLFWLSMDYVASTASI FSVFI LCIDRYRSVQQPLRYLRYRTKTRASATILGAWFFSFLWVIPILGWHHFMPPAPE LREDKC ETDFYNVTWFKIMTAIINFYLPDLLMLWFYVKIYKAVRRHCQHRQLTNGSLPSF SELKLR SDDTKEGAKKPGRESPWGVVKRPSRDPSVGLDQKSTSEDPKMTSPTVFSQE GERETRPCF RLDIMQKQSVAEQDVRGSKANDQALSQPKMDEQSLNTCRRISSETSEDQTLVD QQSFSRTT DSDTSIEPGPGRVKSRSGSNSGLDYIKITWKRLRSHSRQYVSGLHLNRERKAA KQLGFIM AAFILCWIPYFIFFMVI AFCKSCCSEPMHMFTIWLGYINSTLNPLIYPLCNENFKK TFKK ILHIRS</p>
Source	in vitro E.coli expression system
Target Names	Hrh1
Protein Names	Recommended name: Histamine H1 receptor Short name= H1R Short name= HH1R
Expression Region	1-486
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.