



Recombinant Pongo pygmaeus Histamine H1 receptor (HRH1)

Product Code	CSB-CF863085EXP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9N2B0
Product Type	Transmembrane Protein
Immunogen Species	Pongo pygmaeus (Bornean orangutan)
Sequence	<p>MSLPNSSCLLEDKMCEGNKTTMASPQLMPLVVVLSTISLVTVGLNLLVLYAVRS ERKLHT VGNLYIVSLSVADLIVGAVVMPMNILYLLMSKWSLGRPLCLFWLSMDYVASTAS IFSIFI LCIDRYRSVQQPLRYLKYRTKTRASATILGAWFLSFLWVIPILGWNHFRQQISV RREDKC ETDFYDVTWFKVMTAINFYLPDLLMLWFYAKIYKAVQKHCQHRELINGSLPSF SEIKLR PENPKGDAKKPGKESPWEVLKRKPKDAGGGSVLKSPSQTPKEMKSPVFSQ EDDGEVDKL HCFPLDIVQMQTVAEGSSRDYVAINQSHGQLKTDEQGLNTHGASEISEDQML GDSQSFSR TDSDTTETAPGKGLRSGSNTGLDYIKFTWKRLRSHSRQYVSGLHMNRERK AAKQLGFI MAAFILCWIPYFIFFMVIAFCKNCCNEHLHMFTIWLGYINSTLNPLIYPLCNENFK KTFK RILHIRS</p>
Source	in vitro E.coli expression system
Target Names	HRH1
Protein Names	Recommended name: Histamine H1 receptor Short name= H1R Short name= HH1R
Expression Region	1-487
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.