



# Recombinant Pan troglodytes Histamine H1 receptor (HRH1)

<b>Product Code</b>	CSB-CF885636EQV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q9N2B2
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Pan troglodytes (Chimpanzee)
<b>Sequence</b>	<p>MSLPNSSCLLEDKMCEGNKTTMASPQLMPLVVVLSTICLVTVGLNLLVLYAVR  SERKLHT  VGNLYIVSLSVADLIVGAVVMPMNILYLLMSKWSLGRPLCLFWLSMDYVASTAS  IFSIFI  LCIDRYRSVQQPLRYLKYRTKTRASATILGAWFLSFLWVIPILGWNHFMQQTSV  RREDKC  ETDFYDVTWFKVMTAINFYLPDLLMLWFYAKIYKAVRQHCQHRELINGSLSFS  SEIKLR  PENPKGDAKKPGKESPWEVLKRKPKDAGGGSVLKSPSQTXKEMKSPVFSQ  EDDREVDKL  HCFPLDIVPMQTAAEGSSRDYVAVNQSHGQLKTDEQGLNTHGASEISEDQML  GDSQSFSR  TDSDTTETAPGKGKLRSGSNTGLDYIKFTWKRLRSHSRQYVSGLHMNRERK  AAKQLGFI  MAAFILCWIPYFIFFMVIAFCKNCCNEHLHMFTIWLGYINSTLNPLIYPLCNENFK  KTFK RILHIRS</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	HRH1
<b>Protein Names</b>	Recommended name: Histamine H1 receptor Short name= H1R Short name= HH1R
<b>Expression Region</b>	1-487
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>