



Recombinant Drosophila melanogaster Probable G-protein coupled receptor Mth-like 11 (mthl11)

Product Code	CSB-CF306606DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P83118
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	RDIPNCDFFDTVQLRESEKLCNGSYRYEDVVIPAKLTGKYDYEIDYDGDRVSV PKHIRGC VCKLKTCIRFCCHHKLMAGNLCSQDVYENLTYEYTLTDITQLNGSVIKKHVLND MVVQQD LPLPCERHYSLDAETSTYDMWSLYENGLFRHFDQRYLSKQEFCLQPNPTST GKNYSLIV AFNCIQKPSMKMAYGRFECVRKSRLSNASIPVKFSSVFFMVITIAAYLWLPKFR SLHGKC CNLYFICLAITFLNVISLFGIFELKTPICYLTGYAGYFTVMATFLWLSVISFDVWR RFA MRKFQVFYKNKRSSFFNYNIIVWSSAGLLTCIIFLVDQFVETNLDNPYNPAVGV FSCWIF TNGWSATFYFYAPLAILIILNCASFLLTTRYIYVENKQNQKVLNNSSEPQKLSRNH ANYRI YFRLFIIMGGSWFLEIIAFICEMENMWKPLIILNDYINCSQGIIFVATFCNHEMFR LIR KRWVFCDSTP
Source	in vitro E.coli expression system
Target Names	mthl11
Protein Names	Recommended name: Probable G-protein coupled receptor Mth-like 11 Alternative name(s): Protein methuselah-like 11
Expression Region	21-510
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.