



Recombinant Mouse Frizzled-4 (Fzd4)

Product Code	CSB-CF733770MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q61088
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	FGDEEERRCDPIRIAMCQNLGYNVTKMPNLVGHELQTD AELQLTTFTPLIQYG CSSQLQF FLCSVYVPMCTEKINIPIGPCGGMCLSVKRRCEPVLREFGFAWPDTLNCSKFP PQNDHNNH MCMEGPGDEEVPLPHKTPIQPGEECHSVGSNSDQYIWVKRSLNCVLKCGYDA GLYSRSAK EFTDIWMAVWASLCFISTTFTVLTFLIDSSRFSYPERPIIFLSMCYNIYSIAYIVRL TVG RERISCDFEAAEPVLIQEGLKNTGCAIIFLLMYFFGMASSIWWVILTLTWFLAA GLKWG HEAIEMHSSYFHIAAWAIPAVKTIVILIMRLVDADEL TGLCYVGNQNLDALTGFV VAPLF TYLVIGTLFIAAGLVALFKIRSNLQKDGTGKTDKLERLMVKIGVFSVLYTVPATCVI ACYF YEISNWALFRYSADDSNMAVEMLKIFMSLLVGITSGMWIWSAKTLHTWQKCSN RLVNSGK VKREKRGNGWVKPGKGNETVV
Source	in vitro E.coli expression system
Target Names	Fzd4
Protein Names	Recommended name: Frizzled-4 Short name= Fz-4 Short name= mFz4 Alternative name(s): CD_antigen= CD344
Expression Region	37-537
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.