



Recombinant Chicken Frizzled-4 (FZD4)

Product Code	CSB-CF862309CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9IA05
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	FGDEEERRCDAIRIAMCQNLGYNVTKMPNLVGHELQADAELQLTTFTPLIQYG CSSQLQF FLCSVYVPMCTEKINIPIGPCGGMCLSVKRRCEPVLKEFGFAWPDSLNCCKFP PQNDHNNH MCMEGPGDEEVPLHSKTSLQPGEECHSMGSNSDQYIWVKRNLDCVLKCGYD AGLYSRSAK EFTDIWMAVWASLCFISTAFTVLTFLIDSSRFSYPERPIIFLSMCYNIYSIAYIVRL TVG RERISCDFEAAEPVLIQEGLKNTGCAIIFLLMYFFGMASSIWWVILTLTWFLAA GLKWG HEAIEMHSSYFHIAAWAIPAVKTIVILIMRLVDADELTLGLCYVGNQNLDAITGFV VAPLF TYLVIGTLFIAAGLVALFKIRSNLQKDGTGKTKLERLMVKIGVFSVLYTVPATCVI ACYF YEISNWAVFRYSADDSNMAVEMLKIFMSLLVGITSGMWIWSAKTLHTWQKCS NRLVNSGK VKREKRADGWVKPGKGNETVV
Source	in vitro E.coli expression system
Target Names	FZD4
Protein Names	Recommended name: Frizzled-4 Short name= Fz-4 Short name= cFz-4
Expression Region	25-525
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.