



Recombinant Rat Frizzled-4 (Fzd4)

Product Code	CSB-CF865610RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9QZH0
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	FGDEEERRCDPIRIAMCQNLGYNVTKMPNLVGHELQTD AELQLTTFTPLIQYG CSSQLQF FLCSVYVPMCTEKINIPGICGGMCLSVKRRCEPVLKEFGFAWPDSLNC SKFP PQNDHNN MCMEGPGDEEVPLPHKTPIQPGEECHSVGTNSDQYIWKRS LNCVLKCGYDA GLYSRSAK EFTDIWMAVWASLCFISTTFTVLTFIDSSRFSYPERPIIFLSMCYNIYSIAYIVRL TVG RERISCDFEAAEPVLIQEGLKNTGCAIIFLLMYFFGMASSIWWVILTLTWFLAA GLKWG HEAIEMHSSYFHIAAWAIPAVKTIVILIMRLVDADEL TGLCYVGNQSLDALTGFV VAPLF TYLVIGTLFIAAGLVALFKIRSNLQKDGTGKTDKLERLMVKIGVFSVLYTVPATCVI ACYF YEISNWALFRYSADDSNMAVEMLKIFMSLLVGITSGMWWIWSAKTLHTWQKCSN RLVNSGK VKREKRGNGWVKPGKGNETVV
Source	in vitro E.coli expression system
Target Names	Fzd4
Protein Names	Recommended name: Frizzled-4 Short name= Fz-4 Short name= rFz4 Alternative name(s): CD_antigen= CD344
Expression Region	38-538
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.