



Recombinant Arabidopsis thaliana MLO-like protein 3 (MLO3)

Product Code	CSB-CF822027DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q94KB9
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	<p>MTDKEESNHSSEVGAVRSLQETPTWALATVCCFFIAVSICLERLINLLSTRLLKKN RKTSL LEAVEKLSVLMVLGFMSLMLNVTEGEVSKICIKYANRMLPCRKTIKSHNDVSE EDDDD DDGDNHDNSFFHQCSSKGTSLISEEGLTQLSYFFVFLACMHILCNLAILLGM AKMRKW NSWEKETQTVEYLAANDPNRFRITRDTTFARRHLSSWTETSFQLWIKCFFRQF YNSVAKV DYLTLRHGFIFAHVSSNNAFNFQNYIQRSLHEDFKTVVGISPLMWLTVVIFMLLD VSGWR VYFYMSFVPLIIVLVIGTKLEMIVAKMAVTIKENNSVIRGTPLVESNDTHFWFSNP RFL SILHYTLFLNTFEMAFIVWITWQFGINSCYHDNQGIITRLVLAVTVQFLSSYITLP LYA IVTQMGSSYKRAILEEQLANVLRHWQGMVRDKKKTIQTPDNDNNSNNNGDID SGESPVQ TEVASEFRFSGRQSPILQEIQEKTER</p>
Source	in vitro E.coli expression system
Target Names	MLO3
Protein Names	Recommended name: MLO-like protein 3 Short name= AtMlo3
Expression Region	1-508
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.