



Recombinant Mouse Probable G-protein coupled receptor 176 (Gpr176)

Product Code	CSB-CF770482MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q80WT4
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MGHNSSWVSPNTSHPRNTSGAEAGANLSAFGELSEAQLYRQFTTTVQVVIFIG SLLGNFM VLWSTCRTTVFKSVTNRFIKLNACSGICASVVCVPFDIILSSSPHCCWWIYTMLF CKVLK FLHKVFCSVTVLSFPAIALDRYYSVLYPLERKISDAKSRELVMIWAHAVVASVP VFAVT NVADIYAMSTCTEVWSNSLGHLYVLIYNVTTVIVPVAVVFLFLILIRRALASQK KKVI IAALRTPQNTISIPYASQREAEHLHATLLSMVTVFILCSVYATLVVYQTVLNVNPT SVFL LLTAIWLPKVSLLANPVLFLT VNKSVRKCLVGTLVQLHHRYSRRNVVSTGSGVA EPSLEP SMRSGSQLLEMFHIGQQQIFKPSEDEEESEAKYLGSAADFQAKEVLTSCPEGE QEPPQLAP SVPPPGTVDSEPRVSPVAPMESGIFPKYSLQFGFGPFELPPQWLSETRNSK KRLLPPLG NTPEELIQTKVPRVNRVERKMSRNNKVSIFPKVDS
Source	in vitro E.coli expression system
Target Names	Gpr176
Protein Names	Recommended name: Probable G-protein coupled receptor 176 Alternative name(s): G-protein coupled receptor AGR9
Expression Region	1-515
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.