



Recombinant Chicken Frizzled-2 (FZD2)

Product Code	CSB-CF884854CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9IA06
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	<p>PDHGFCQPISIPLECTDIAYNQTIMPNLLGHTNQEDAGLEVHQFYPLVKVQCSLE LKFFLC SMYAPVCTVLEQAIPPCRSICERARQGCEALMNKFGFQWPERLRCENFPRHG AEQICVGQ NHSEGGSSAALLTSAAPPAAHGTPGAPRYATPDRPFHCPRALKVPGYLYNK FLGEKDCA APCEPSRPDGHMFFNEDEIRFARVWILVWSVLCCASTFFTVTTYLVDMQRFY PERPIIF LSGCYTMVSVAYIAGFVLEERVVCNERFQEDGYRTVVQGTKKEGCTILFMMMLY FFSMAS IWWWILSLTWFLAAGMKWGHEAIEANSQYFHAAWAVPAVKTITILAMGQIDGD LLSGVC FVGLNGIDPLRGFVLAPLFVYLFIGTSFLLAGFVSLFRIRTIMKHGGTKTEKLERL MVRI GVFSVLYTVPATIVIACYFYEQAFRQHWERSWISQHCPSLAIPCPLHFTPRMTP DFTVYM IKYLMTLIVGITSGFWIFSGKTLHSWRKFYTRLNSRQGETTV</p>
Source	in vitro E.coli expression system
Target Names	FZD2
Protein Names	Recommended name: Frizzled-2 Short name= Fz-2 Short name= cFz-2
Expression Region	1-523
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.