



Recombinant *Drosophila sechellia* Protein wntless (wls)

| | |
|--------------------------|---|
| Product Code | CSB-CF451838DMG |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | B4HEB1 |
| Product Type | Transmembrane Protein |
| Immunogen Species | <i>Drosophila sechellia</i> (Fruit fly) |
| Sequence | APLPAGHVTVLGSLCREDHARQNDTGFLLYSRGAGACIPVTREEVEQDSTKM ANELVHVF QMPLPRDLRDLDYSRWQQNLIGVLQVEFGYDSSSELREPPRELQLTIDMRLAY RNKGDPD NGWKLYAHGVEHRYLDCVTSHVGPTEPLYSCDMIPLFELGALHHSFYLLNLR PLDTPSQ MNLQFGHMHDLTLTAIHQNGGFTQIWLLKTMFLPFVVGIMIWFWRRVHLLQR SPALLEY MLIYLGAAALTFLNLPLEYLSLVYEMPYMLLLSDIRQGIFYAMLLTFWLVFAGEHM LIQDA PNKSTIRSRYWKHLSAVVVGICISLFVFDICERGVQLRNPFYSIWTTPLGAKVAM TFIVLA GVSAAIYFLFLCYMIWKVFRNIGDKRTSLPSMSQARRLHYEGLIYRFKFLMLATL VCAAL TVAGFIMGQMAEGQWDWNDNVAIQPTSAFLTGVYGMWNIYIFALLILYAPSHK QWPAMHH SDETTQSNENIVASAASEEIEFSHLPSDSNPSEISSLTSFTRKVAFD |
| Source | in vitro E.coli expression system |
| Target Names | wls |
| Protein Names | Recommended name: Protein wntless |
| Expression Region | 36-562 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |