



Recombinant *Drosophila erecta* Protein wntless (wls)

Product Code	CSB-CF459258DLK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3NGS7
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila erecta</i> (Fruit fly)
Sequence	APLPAGHVTVLGSLCREDHGRQNDTGFLLYSRGAGACIPVTREEVEQDSTKM ANELVHVF QMPLPRDLRDLDYSRWQQNLIGVLQVEFGYDSSSELREPPRELQLTIDMRLAY RNKGDPD NGWKLYAHGVEHRYLDCVSSHVGPTEETLYSCDMIPLFELGALHHSFYLLNLRP PLDTPSQ MNLQFGHMHDLTLTAIHQNGGFTQIWLLKTVLFPFVVGIMIWFWRVHLLQR SPALLEY MLIYLGAAALTFLNLPLEYLSLVFEMPYMLLLSDIRQGIFYAMLLTFWLVFAGEHM LIQDA PNKSTIRSRYWKHLSAVVVGICISLFFVICERGVQLRNPFIYSIWAMPLAAKMA MTFIVLA GVSAAIYFLFLCYMIWKVFRNIGDKRTSLPSMSQARRLHYESLIYRFKFLMLATI VCAAL TVTGFIMGQRAEGQWDWNDNVAIQPTSAFLTGVYGMWNIYIFALLILYAPSHK QWPTMHH SDETTQSNENIVASAASEEIEFSHLPSDSNPSEISSLTSFTRKVAFD
Source	in vitro E.coli expression system
Target Names	wls
Protein Names	Recommended name: Protein wntless
Expression Region	36-562
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.