



Recombinant Drosophila melanogaster Putative gustatory receptor 66a (Gr66a)

Product Code	CSB-CF895186DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VSH2
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MAQAEDAVQPLLQQFQQLFFISKIAGILPQDLEKFRSRNLEKSRNGMIYMLST LILYVV LYNILIYSFGEDRSLKASQSTLTFVIGLFLTYIGLIMMVSDQLTALRNQGRIGEL YERI RLVDERLYKEGCVMDNSTIGRRIRIMLIMTVIFELSILVSTYVKLVDYSQWMSLL WIVSA IPTFINTLDKIWFVAVSLYALKERFEAINATLEELVDTHEKHKLWLRGNQEVPPPL DSSQP PQYDSNLEYLYKELGGMDIGSIGKSSVSGSGKNKVAPVAHSMNSFGEAIDAAS RKPPPPP LATNMVHESELGNAAKVEEKLNNLCQVHDEICEIGKALNELWSYPILSLMAYGF LIFTAQ LYFLYCATQYQSIPSLFRSAKNPFITVIVLSYTSKGKCVYLIYLSWKTSQASKRTGI SLHK CGVVADDNLLYEIVNHLNLSKLLNHSVDFSACGFFTLDMETLYGVSGGITSYLILIL QFNL AAQQAKEAIQTFNSLNDTAGLVGAATDMDNISSTLRDFVTTTMTPAV
Source	in vitro E.coli expression system
Target Names	Gr66a
Protein Names	Recommended name: Putative gustatory receptor 66a
Expression Region	1-527
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.