



Recombinant Bovine G protein-coupled receptor 161 (GPR161)

Product Code	CSB-CF009774BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2YDN1
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MSLNSSLGHRKELSNLTEGASDQGGSGVTEFVAIVIITVFVCLGNLVIVITLYRK SYLLT LSNKFVFSLTLSNFLLSVLVLPFVVTSSIRREWIFGVVWCNFSALLYLLISSASML TLGI IAVDRIYYAVLYPMAYPMKITGNRAVMVLAYIWLHSLIGCLPPLFGWSSVEFDEF KWMCA AWHREPGYTAFWQIWCALLPFLVMLVCYGFIFRVARVKARKVHCGAVVTVEV GVQRTGRK NSSTSTSSSGSRKSAFQGVVYSANQCKALVTILVVIGAFMVTWGPYMVVITSE ALWGKNC VSPTLETWATWLSFTSAICHPLIYGLWNKTVRKELLMCFGDRYYREPFVQRQ RTSRLFS ISNRITDLGLSPHLTALMAGEQPLGNSSSTGDTGFSCSQDSGTDVMLLEDYTS DDNPLHG TCPPKRRSSVTFEDEVETQIKEAAKNPILHVKADVHKSLSYATSLAKAIEAEAKI NLFGE EALPGVLLTARTVPGIGFGSRRGRSRTLAGQRLQLQSIEEGDVLATEQR
Source	in vitro E.coli expression system
Target Names	GPR161
Protein Names	Recommended name: G protein-coupled receptor 161
Expression Region	1-528
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.