



Recombinant Mouse Frizzled-7 (Fzd7)

Product Code	CSB-CF730721MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q61090
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	QPYHGEKGISVPDHGFCQPISIP LCTDIAYNQ TILPNLLGHTNQEDAGLEVHQF YPLVKV QCSPELRFFLC SMYAPVCTVLDQAIPPCRSLCERARQGCEALMNKFGFQWPE RLRCENFP VHGAGEICV GQNTSDGSGGAGGSPTAYPTAPYLPDPPFTAMSPSDGRGRLSF PFSCPRQL KVPPYLG YRFLGERDCGAPCEPGRANGLMYFKEEERRFARLWVGWVSVLSC ASTLFTVLT YLVDMMRRFSYPERPIIFLSGCYFMVAVAHVAGF LLEDRAVCVERFSDDGYRTV AQGTKKE GCTILFMVLYFFGMASSIWWVILSLTWFLAAGMKWGHEAIEANSQYFH LAAWA VPAVKTI TILAMGQVDGDLLSGVCYVGLSSVDALRGFVLAPLFVYLFIGTSFLLAGFVSLF RIRTIM KHDGTEKLEKLMVRIGVFSVLYTVPATIVLACYFYEQAFREHWERTWLLQT CKSYAVP CPPRHFS PMSPDFTVFMIKYLMTMIVGITTGFWIWSGKTLQSWRRFYHRLSHS SKGETAV
Source	in vitro E.coli expression system
Target Names	Fzd7
Protein Names	Recommended name: Frizzled-7 Short name= Fz-7 Short name= mFz7
Expression Region	33-572
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.