



Recombinant Mouse Frizzled-2 (Fzd2)

Product Code	CSB-CF881309MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9JIP6
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	QFHGEKGISIPDHGFCQPISIP LCTDIAYNQTIMPNLLGHTNQEDAGLEVHQFYPLVKVQ CSPELRFFLC SMYAPVCTVLEQAIPPCR SICERARQGCEALMNKFGFQWPERLRCEHFPR HGAEQICV GQNHSE DGAPALLTTAPPSGLQPGAGGTPGGPGGGGSPPRYATLEHPPHCPR VLKVPSYLSYKFLGERDCAAPCEPARPDGSMFFSQEETRFARLWILTWSVLCCASTFFTV TTYLVDMQRFRYPERPIIFLSGCYTMVSVAYIAGFVLQERVVCNERFSEDGYRTVVQGTK KEGCTILFMMLYFFSMASSIWWVILSLTWFLAAGMKWGHEAIEANSQYFHLLAAWAVPAVK TITILAMGQIDGDLLSGVCFVGLNSLDPLRGFVLAPLFVYLFIGTSFLLAGFVSLFRIRT IMKHDGKTEKLERLMVRIGVFSVLYTVPATIVIACYFYEQAFREHWERSWVSQHCKSLA IPCPAHYTPRMSPDFTVYMIKYLMTLIVGITSGFWIWSGKTLHSWRKFYTRLTNSRHGET TV
Source	in vitro E.coli expression system
Target Names	Fzd2
Protein Names	Recommended name: Frizzled-2 Short name= Fz-2 Short name= mFz2 Alternative name(s): Frizzled-10 Short name= Fz-10 Short name= mFz10
Expression Region	29-570
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.