



Recombinant Human Frizzled-7 (FZD7)

Product Code	CSB-CF009110HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75084
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>QPYHGEGKISVPDHGFCQPISIP LCTDIAYNQ TILPNLLGHTNQEDAGLEVHQF YPLVKV QCSPELRFFLC SMYAPVCTVLDQAIPPCRSLCERARQGCEALMNKFGFQWPE RLRCENFP VHGAGEICV GQNTSDGSGGGPGGGPTAYPTAPYLPDLPFTALPPGASDGRGRP AFPFSCPR QLKVPPYLG YRFLGERDCGAPCEPGRANGLMYFKEEERRFARLWVGWWSVL CCASTLFTV LTYLVD MRRFSYPERPIIFLSGCYFMVAVAHVAGF LLEDRAVCVERFSDDGYR TVAQGTK KEGCTILFMVLYFFGMASSIWWVILSLTWFLAAGMKWGHEAIEANSQYFH LAA WAVPAVK TITILAMGQVDGDL LSGVCYVGLSSVDALRGFVLAPLFVYLFIGTSFLLAGFVSL FRIRT IMKHDG TKTEKLEKLMVRIGVFSVLYTVPATIVLACYFYEQAFREHWERTWLLQ TCKSYA VPCPPGHFPPMSPDFTVFM IKYLMTMIVGITTGFWIWSGKTLQSWRRFYHRLS HSSKGET AV</p>
Source	in vitro E.coli expression system
Target Names	FZD7
Protein Names	Recommended name: Frizzled-7 Short name= Fz-7 Short name= hFz7 Alternative name(s): FzE3
Expression Region	33-574
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.