



Recombinant Mouse G-protein coupled receptor 161 (Gpr161)

Product Code	CSB-CF009774MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B2RPY5
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p>MDFVQHALLTASRGALTMSLNSSLSYRKELSNLTATEGGEGGAVSEFIAIIITVL VCLG NLVIVVTLYKKSYPYLLTSLNKFVFSLTLSNLLSVLVLPFVVTSSIRREWIFGVVWC NFSA LLYLLISSASMLTLGVIAIDRYAVLYPMVYPMKITGNRAVMALVYIWLHSLIGCL PPLF GWSSVEFDEFKWMCVAAWHQEPGYTIFWQIWCALFPFLIMLVCYGFIFRVAR VKARKVHC GTVVTVEEDSQRSGRKNSSSTSTSSSGSRRNALQGVVYSANQCKALITILVVIG AFMVTWG PYMVVITSEALWGKNCVSPTLETWATWLSFTSAICHPLIYGLWNKTVRKELLG MCFGDRY YRESFVQRQRTSRLFSISNRITDLGLSPHLTALMAGGQSLGHSSSTGDTGFSY SQDSGTD VMLLEDGTSEDNPPQHCTCPPKRRSSVTFEDEVQIKEAAKNSLLHVKAIEVHK SLDSYAA SLAKAIEAEAKINLFGEEALPGVLFTARTVPGAGFGGRRGSRTLNVNQLQLQSI KEGNVLAAEQR</p>
Source	in vitro E.coli expression system
Target Names	Gpr161
Protein Names	Recommended name: G-protein coupled receptor 161
Expression Region	1-545
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.