



Recombinant Calcitonin receptor-like protein 1 (pdf-1)

Product Code	CSB-CF600242CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q09460
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MADATSPFNVSILDNSTKLSEMVESGWNVLASTSVQAFNEAMDVLEESYPLCK KMLDHNN LFPERDPNDTRIWCNATYDTVLCWPPTPANSSVTLQCPHMKGLDPNKNITKD CHVSGVWS GRNAGEMGPTLPGWTNFTMCYTDEVIYIMQNLNNESTIAQEVARNARKLEFV GLGLSLV SLILAISIFSFRRLRVFRNLLHLHLMIAMLMVVILRLVLYIDLIFTGENGPHTNSA EGK TINTMPIVCEGMFFFLEYFKTVTFCWMFLEGIYLNQIVFGFFNSEPKLLPYFIA GYGIP LVHTMLWLLVVLIKKDFKVERCLGSYYLEPEFWILDGPRMAELVINLFFICNVIR VLYSK VRESNNTSEAGLKKS VKAAMMLLPLLGPVNIMQTIPFAPTRDNIMVFAVWTYTA SFTYMY QGLMVASIYCFTNKEVNHVLKTFYARYRLLHKSQNELRRGSRVASHYAAKN GTANASAP QTNNADEFGKLSPPFSRSKKGSDSTTKLMKDAVMEEKKNANNNGYGSAGE MTPLREGSN RSTKSP
Source	in vitro E.coli expression system
Target Names	pdf-1
Protein Names	Recommended name: Calcitonin receptor-like protein 1 Alternative name(s): PDF Receptor homolog 1
Expression Region	1-546
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.