



# Recombinant Escherichia coli Thiol:disulfide interchange protein DsbD (dsbD)

<b>Product Code</b>	CSB-CF633364EGW
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q1R3C4
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Escherichia coli (strain UTI89 / UPEC)
<b>Sequence</b>	<p>GLFDAPGRSQFVPADQAFDFQQNQHDLNLTWQIKDGYLYRKQIRITPEHA            KIADVQL            PQGVWHEDEFYKGSEIYRDRLTLPVTINQASAGATLTVTYQGCADAGFCYPPE            TKTVPLS            EVVANNEASQPVSVPPQQEQPTAQLPFSALWALLIGIGIAFTPCVLPMPYPLISGIV            LGGKQ            RLSTARALLLTFIYVQGMALTYTALGLVVAAGLQFQAALQHPYVLIGLAIVFTLL            AMSM            FGLFTLQLPSSLQTRLTMSNRQQGGSPGGVFIMGAIAGLICSPCTTAPLSAILL            YIAQS            GNMWLGGGTLYLYALGMGLPLMLITVFGNRLLPKSGPWMEQVKTAFGFVILAL            PVFLLER            VIGDIWGLRLWSALGVAFFGWAFITSLQAKRGWMRVVQIILLAAALVSVRPLQD            WAFGET            HTAQTQTHLNFTQIKTVDELNQUALVEAKGKPVMLDLYADWCVACKEFEKYTFS            DPQVQKA            LADTVLLQANVTANDAQDVALLKHLNVLGLPTILFFDGQQQEHPQARVTGFMD            AETFSAH LRDRQP</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	dsbD
<b>Protein Names</b>	Recommended name: Thiol:disulfide interchange protein DsbD EC= 1.8.1.8 Alternative name(s): Protein-disulfide reductase Short name= Disulfide reductase
<b>Expression Region</b>	20-565
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.