



Recombinant Rat Glucagon-like peptide 2 receptor (Glp2r)

Product Code	CSB-CF009515RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9Z0W0
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MRPQPSPAVPSRCREAPVPRVRAQPVGIPEAQGPVPLHSQQMRLWGPGRP FLALLLLVS IKQVTGSELLKETTQKWANYKEKCLEDLHNRLSGIFCNGTFDRYVCWPHSYPG NVSVPCPS YLPWWNAESPGRAYRHCLAQGTWQTRENTTDIWWQDESECSSENHSFRQNVD HYALLYTLQL MYTVGYSVSLISLFLALTLFLFLRKLHCTRNYIHMNLFASFILKVLAVLVKDMVSH NSYS KRPDESSEGWMSYLSETSVSCRSVQVLLHYFVGTNHLWLLVEGLYLHTLLEPT VFPERRLW PKYLVVVGWAFPMFLVIPWGFARAHLENTRCWATNGNLKIWWIIRGPMLLCVTV NFFIFLK ILKLLISKLKAHQMCFRDYKYRLAKSTLLLIPLLVGHEVLFTFFPDDQVQGFSKRI RLF QLTLSSVHGFLVALQYGFANGEVKAELRKSWSGRFLLARHWGCRTCVLGKNFR FLGKCSKK LSEGDSSETLQKLRFFSTCSSHLASETLGDVGVQPHRGRGAWPRGSSLSSESSE GDFTLANT MEEILEESEI
Source	in vitro E.coli expression system
Target Names	Glp2r
Protein Names	Recommended name: Glucagon-like peptide 2 receptor Short name= GLP-2 receptor Short name= GLP-2-R Short name= GLP-2R
Expression Region	1-550
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.