



# Recombinant *Drosophila melanogaster* Probable G-protein coupled receptor Mth-like 10 (mthl10)

<b>Product Code</b>	CSB-CF893356DLU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q9W0R5
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Sequence</b>	<p>IPFGTYVMAERDHYHTIDDPNVPCNFYDTVNLTGHRLFNGSYDYYGTIVPAEL            VGTYDY            IHSSLTERIEVREHVRGCVCKFKSCLNICCPWRQVFNSEVDGCIIDHSDNRTW            PDPPMLN            IFRNESTILVNMFTQFAIQSFRPCPKMFSLQPETNNWDDYLLFENGSMRLRVD            DKLLIRK            NEFCMVPTYVNESDMFYTIHPANCDMQDDHSTVKIINSYAMMFSIPFMMLTIAV            YLLIPE            LRNQHGKSLVCYLIGLSVGYSSLCYVQLYQVDATGVTCKVFGYTAYFFFMGAY            MWLSVIS            FDLWHNFRGTRGINRFQEKKRFLFYSLYSWGIALVFLAFTYCAQQLTNLPANL            KPGIGDG            VYCWLDMSNWAAMIYFYGPILAIVVANTIMFIMTAIKIHGVQREMARIASENSTK            NLRT            EKDKRFYRAWSNYRFGFLRLFLIMGITWLTELISYFVGSDKGWSKLFYISDLA            NAMQGF            LIFMLFVMKKKVKHLITNRCSSVRDGSNQRQSQYSTKTTSSSVANLSLHEKPS            VEKPLVI SSSVDPQKTTIFR</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	mthl10
<b>Protein Names</b>	Recommended name: Probable G-protein coupled receptor Mth-like 10 Alternative name(s): Protein methuselah-like 10
<b>Expression Region</b>	33-585
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.