



Recombinant *Xenopus laevis* Frizzled-8 (fzd8)

Product Code	CSB-CF009111XBE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O93274
Product Type	Transmembrane Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Sequence	<pre>KELSCQEITVPLCKDIGYNYTYMPNQFNHDTQDEAGMEVHQFWPLVVIHCSPD LKFFLCS MYTPICLEDYKKPLPPCRSVCERARAGCAPLMRQYGFAPDRMRCDRLEPEQ GNPDTLCMD YYNRTEQTTAAPSHPEPPKPPARSVPKGRTRVEPPRSRSRATGCESGCQCR APMVQVSNE RHPLYNRVRTGQIPNCAMPCHNPFFSPEERTFTEFWIGLWSVLCFASTFATVS TFLIDME RFKYPERPIIFLSACYLLVSTGYLIRLIAGHEKVACSRGELDLEHIIHYETTGPALC TLV FLLIYFFGMASSIWWVILSLTWFLAAGMKWGNEAIAGYSQYFHAAWLVPSTIKS IAVLAL SSVDGDPVAGICFVGNQNLNLRGFVLAPLVIYLFIGSMFLLAGFVSLFRIRSVI KQGGT KTDKLEKLMIRIGIFSVLYTVPATIVVACFFYEQHNRRQGWEVAHNCNSCQPEMA QPHRPD YAVFMLKYFMCLVVGITSGVWIWSGKTLSEWRAFCTRCCWGSKATGGSMYS DVSTGLTWR SGTGSSVSCPQMPLSQV</pre>
Source	in vitro E.coli expression system
Target Names	fzd8
Protein Names	Recommended name: Frizzled-8 Short name= Fz-8 Short name= Xfz8
Expression Region	24-581
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.