



Recombinant Human Frizzled-5 (FZD5)

Product Code	CSB-CF613403HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q13467
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>ASKAPVCQEITVPMCRGIGYNLTHMPNQFNHDTQDEAGLEVHQFWPLVEIQCS SPDLRFFL CSMYTPICLPDYHKPLPPCRSVCERAKAGCSPLMRQYGFAPWPERMSCDRLP VLGRDAEVL CMDYNRSEATTAPPRPFPKPTLPGPPGAPASGGECPPAGGPFVCKCREPFVP ILKESHPL YNKVRTGQVPNCAVPCYQPSFSADERTFATFWIGLWSVLCFISTSTTVATFLID MERFRY PERPIIFLSACYLCVSLGFLVRLVVGHASVACSREHNHIHYETTGPALCTIVFLLY YFFG MASSIWWVILSLTWFLAAGMKWGNEAIAQYFHLAAWLIPSVKSITALALSS VDGPV AGICYVGNQNLNSLRGFVLGPLVLYLLVGTFLLAGFVSLFRIRSVIKQGGTKTD KLEKL MIRIGIFLLYTPASIVVACYLYEQHYRESWEAALTCACPGHDTGQPRAKPEY WVLMLK YFMCLVVGITSGVWIWSGKTVESWRRFTSRCCCRPRRGHKSGGAMAAGDYP EASAALTGR TGPPGPAATYHKQVSLSHV</p>
Source	in vitro E.coli expression system
Target Names	FZD5
Protein Names	Recommended name: Frizzled-5 Short name= Fz-5 Short name= hFz5 Alternative name(s): FzE5
Expression Region	27-585
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.