



Recombinant Rat Frizzled-5 (Fzd5)

Product Code	CSB-CF807066RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8CHL0
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>ASKAPVCQEITVPMCRGIGYNLTHMPNQFNHDTQDEAGLEVHQFWPLVEIHC SPDLRFFL CSMYTPICLPDYHKPLPPCRSVCERAKAGCSPLMRQYGFAPWPERMSCDRLP VLGDAEVL CMDYNRSEATTASPKSFPKPTLPGPPGAPSSGGECPSGGPSVCTCREPFVP ILKESHPL YNKVRTGQVPNCAVPCYQPSFSPDERTFATFWIGLWSVLCFISTSTTVATFLID MERFRY PERPIIFLSACYLCVSLGFLVRLVVGHASVACSREHSHIHYETTGPALCTVVFL VYFFG MASSIWWVILSLTWFLAAGMKWGNEAIAQYFHLAAWLIPSVKSITALALSS VDGPV AGVCYVGNQNLNSLRGFVLGPLVLYLLVGTFLLAGFVSLFRIRSVIKQGGTKT DKLEKL MIRIGIFLLYTPASIVVACYLYEQHYRESWEAALTCACPGSDAGQPRAKPEY WVLMLK YFMCLVVGITSGVWIWSGKTLESWRRFTSRCCSSRRGHKSGGAMVAGDYA EASAALTGR TGPPGPAAAYHKQVSLSHV</p>
Source	in vitro E.coli expression system
Target Names	Fzd5
Protein Names	Recommended name: Frizzled-5 Short name= Fz-5
Expression Region	27-585
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.