



Recombinant *Drosophila melanogaster* Protein wntless (wls)

Product Code	CSB-CF850199DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q95ST2
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	APLPAGHVTVLGSLCREDHARQNDTSFLLYSRGAGACIPVTREEVEQDSTKM ANELVHVF QMPLPRDLRDLDYSRWQQNLIGVLQVEFGYDSSSELREPPRELQLTIDMRLAY RNKGDPD NGWKLYAHGVEHRYLDCVTSHVGPTEPLYSCDMIPLFELGALHHSFYLLNLR PLDTPSQ MNLQFGHMHDLTLTAIHQNGGFTQIWLLKTMFLPFVVGIMIWFWRRVHLLQR SPALLEY MLIYLGAAALTFLNLPLEYLSLVYEMPYMLLLSDIRQGIFYAMLLTFWLVFAGEHM LIQDA PNKSTIRSRYWKHLSAVVVGICISLFVFDICERGVQLRNPFYSIWTTPLGAKVAM TFIVLA GVSAAIYFLFLCYMIWKVFRNIGDKRTSLPSMSQARRLHYEVPLDQKVEDWAG IVYFYTK AFFFQLHKANESKGLIYRFKFLMLATLVCAALTVAGFIMGQMAEGQWDWNDN VAIQPTSA FLTGVYGMWNIYIFALLILYAPSHKQWPTMHHSDETTQSNENIVASAASEEIEF SHLPSD SNPSEISLTSFTRKVAFD
Source	in vitro E.coli expression system
Target Names	wls
Protein Names	Recommended name: Protein wntless Alternative name(s): Evenness interrupted Sprinter
Expression Region	36-594
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.