



# Recombinant Didelphis marsupialis virginiana Parathyroid hormone/parathyroid hormone-related peptide receptor (PTH1R)

<b>Product Code</b>	CSB-CF018988DKL
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P25107
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Didelphis virginiana (North American opossum) (Didelphis marsupialis virginiana)
<b>Sequence</b>	DADDVITKEEQIILLRNAQAQCEQRLKEVLRVPELAESAKDWMRSRAKTKKEK PAEKLYS QAEESREVSQRSRLQDGFCLPEWDNIVCWPAGVPGKVVAVPCPDYIYDFNHK GRAYRRCD SNGSWELVPGNNRTWANYSECVKFLTNETREREVFDRLGMIYTVGYSISLGS LTVAVLIL GYFRRLHCTRNYIHMHLFVSFMLRAVSIFIKDAVLYSGVSTDEIERITEEELRAF TEPPP ADKAGFVGCRVAVTVFLYFLTNNYYWILVEGLYLHSLIFMAFFSEKKYLWGFTL FGWGLP AVFVAVWVTVRATLANTECWDLSSGNKKWIIQVPILAAIVVNFILFINIIRVLATKL RET NAGRCDTRQQYRLLKSTLVLMPLFGVHYIVFMATPYTEVSGILWQVQMHYE MLFNSFQG FFVAIYCFCNGEVQAEIKKSWSRWTLALDFKRKARSGSSTYSYGPMVSHTSV TNVGPRG GLALSLSPLAPGAGASANGHHQLPGYVKHGSISENSLPSSGPEPGTKDDGY LNGSGLYE PMVGEQPPPLLEEERETVM
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	PTH1R
<b>Protein Names</b>	Recommended name: Parathyroid hormone/parathyroid hormone-related peptide receptor Alternative name(s): PTH/PTHrP type I receptor Short name= PTH/PTHr receptor Parathyroid hormone 1 receptor Short name= PTH1 receptor
<b>Expression Region</b>	27-585
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .