



Recombinant Mouse Frizzled-10 (Fzd10)

| | |
|--------------------------|--|
| Product Code | CSB-CF803892MO |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | Q8BKG4 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Sequence | ISSMDLERPGDGKQCQPVEIPMCKDIGYNTTRMPNLMGHENQREAAIQLHEFAP LVEYGCH SHLRFFLCSLYAPMCTEQVSTPIPACRVMCEQARLKCSPIMEQFKFRWPDSL CSKLPNK NDPNYLCMEAPNNGSDEPSRSGMFPPLFRPQRPHSAQEHLKDGPGRA GCDNPGKFHH VEKSESCAPLCTPGVDVYWSRDDKRFVWVLAWSVLCFFSSAFTVLTFLIDP SRFRYPE RPIIFLSMCYCVYVSVGYIIRLFAGAESIACDRDSGQLYVIQEGLESTGCTLVFLV YYFG MASSLWVVLTLTWFLAAGKKWGHEAIEANSSYFHAAWAIPAVKILILVMR RVAGDEL TGVCYVGSMDVNALTGFVLVPLACYLVIGTSFILSGFVALFHIRVMKTGGENT DKLEKL MVRIGVFSLLYTVPATCVIACYFYERLNMDYWKMLATQHKCKMNNQTKTPDCL MTTSIPA VEVFMVKVSMMLLVGITSGVWVWTSKTLQSWQHVC SRGLKRKSRKSPASVV TSAGIYKKA QHPQKPHLGKYELPAQPSACV |
| Source | in vitro E.coli expression system |
| Target Names | Fzd10 |
| Protein Names | Recommended name: Frizzled-10 Short name= Fz-10 Alternative name(s): CD_antigen= CD350 |
| Expression Region | 22-582 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |