



Recombinant Chicken Frizzled-10 (FZD10)

Product Code	CSB-CF886116CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9PWH2
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	ISSIDIERPGDGRCQPIEIPMCKDIGYNMTRMPNLMGHENQREAAIQLHEFAPL VEYGCH GHLKFFLCSLYAPMCTEQVSTPIACRVMCEQARLKCSPIMEQFNFKWPDSL CSKLPNK NDPNYLCMEAPNNGSDEPPRGSSMLPPMFRPQRPSTGHDLQQHKDSLRTS CENPGKFHH VEKSASCAPLCTPGVDVYWSKDDKQFAVIWIAIWSILCFSSAFTVLTFLIDPQR FKYPE RPIIFLSMCYCVYSVGYIIRLFSGAESIACDRDSGQLYVIQEGLESTGCTIVFLVL YYFG MASSLWWVILTTLTWFLAAGKKWGHEAIEANSSYFHLAAWAIPAVKTIMILVMRR VAGDEL TGLCYVGSMDVNALTGFVLIPLACYLIIGTSFILSGFVALFHIRRVMTGGENTD KLEKL MVRIGVFSVLYTVPATCVIACYFYERLNMDYWKIVASQQKCKMNNQTKNLDC MMNNSIPA VEIFMVKIFMLLVGITSGMWIWTSKTLQSWQNVCSRRLKKRSRRKPASVITSS GIYKKP QHPQKTHLAKYESTLQPPTCV
Source	in vitro E.coli expression system
Target Names	FZD10
Protein Names	Recommended name: Frizzled-10 Short name= Fz-10 Short name= cFz-10
Expression Region	25-585
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.