



Recombinant Human Frizzled-10 (FZD10)

Product Code	CSB-CF891782HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9ULW2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	ISSMDMERPGDGKCKQPIEIPMCKDIGYNMTRMPNLMGHENQREAAIQLHEFAP LVEYGCH GHLRFFLCSLYAPMCTEQVSTPIPACRVMCEQARLKCSPIMEQFNFKWPDSL CRKLPNK NDPNYLCMEAPNNGSDEPTRGSGLPPLFRPQRPHSAQEHPLKDGPPGRGG CDNPGKFHH VEKSASCAPLCTPGVDVYWSREDKRFVWVWLAIVAVLCFFSSAFTVLTFLIDP ARFRYPE RPIIFLSMCYCVYVSVGYLIRLFAGAESIACDRDSGQLYVIQEGLESTGCTLVFLV LYYFG MASSLWVVLTLTWFLAAGKKWGHEAIEANSSYFHAAWAIPAVKTILILVMR RVAGDEL TGVCYVGSMDVNALTGFVLIPLACYLVIGTSFILSGFVALFHIRRVMKTGGENTD KLEKL MVRIGLFSVLYTVPATCVIACYFYERLNMDYWKILAAQHKCKMNNQTKTLDC MAASIPA VEIFMVKIFMLLVGITSGMVIWTSKTLQSWQQVCSRRLKKSRRKPPASVITSG GIYKKA QHPQKTHHGKYEIPAQSPTCV
Source	in vitro E.coli expression system
Target Names	FZD10
Protein Names	Recommended name: Frizzled-10 Short name= Fz-10 Short name= hFz10 Alternative name(s): FzE7 CD_antigen= CD350
Expression Region	21-581
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.