



Recombinant Dictyostelium discoideum Frizzled and smoothened-like protein K (fslK)

Product Code	CSB-CF700832DKK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q54PA5
Product Type	Transmembrane Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Sequence	<p> QQYYPIDPTGKCEQYIGDSVITPCNSIYVSSTANQSTAMLSLNLYFSLGGSSA SCQNAY TFATLCSTYLPECEIFIDNSTNKEIAMPKRVCLDTCNSVSKNCQITSYFDCNQLE PISKL PLFPKESSDYNLTTYGGELNYRLDCYDPPKSNETTIVIGYCPFPLVFRNRTIVGK GSDSLD NSENYGYTYVSEDSYCTIPCPAPIFSEKQWDNIFDTSDAISLVSLLCVYLFITY MVINP KRNKYDYFFSFFVLSIILMSIAGTIGFSVGGTRKLLCPEINRRGVYTDPAVAAAG WIFQF AIINAILWFSINSFELWFQIKFIKRKLHLIKFYILAVLVISIALSVPLSAIGEFNAGLG N FVWIESGKYQNWFFWGPLGIVLTVGTTFIGLVIWEIYKIVSSTNKSDFFKLQLK PLMNM LLIYLTFVYLFGYNFYIHNSLNGFYGSSEEFKNCIISTDGKDCRIQGPPYSSILMF VFCL RIYGVYICIALYGFSPKTRSIWSNSIVFNNSMANKLKTLYISSGTTKGGTSSTDIK MSTNN NSNMDSGGGKSSSMEPDEIILR </p>
Source	in vitro E.coli expression system
Target Names	fslK
Protein Names	Recommended name: Frizzled and smoothened-like protein K
Expression Region	19-580
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>