



Recombinant *Drosophila melanogaster* Frizzled-3 (fz3)

Product Code	CSB-CF009106DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O77438
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	<p>ANGAGHNGPVASGAGPNGLQCQPIAVSACQGLGYNMTALPNLAGHTNQLEA ELQIAKLVP LIESGCSRRARFLLCSSLFPLCTPDVPRPVAACKLLCETVRGECMENAPPELM ELWPSFL NCDGLPQPEKHELQIPQEVAVPGGSPSGPPTTGSPGVEDHPQTYRFWKS GASPTSCLA GVLCPQNFSGSPFNPEECVPQCQRDAFHTSSQKKTSETLILGLSAVCFVLTFLF ALVTFWA EPTRFGYPERPVLFLCLCYNLFSVCYLERIVFHNQARMHDVELQGRLMRPGCL LTPPCLA SYITTSYLSLCAASWWLIFALCFYLSSHKKWSSEALEKRSGLFHVLAWVPPLAP PIAALL LEKVRPSELTGMCYAPGFVELPALVLLLLGLYFTLRASRLLSLQQQLQPTLAH HRFQGI RKRFVLFSLLYFAPTTAGVVAALCERYADSVPCSTPDDCLSPTPLSAWPALV RIFFQLV GGTLTGLWVWSRKTCESYRNRLGASGTPTSSLMNQSKAAGALPKKHLYTSG KSMLPTGGI TPLYAGISFHNVPVYNPNQSRV</p>
Source	in vitro E.coli expression system
Target Names	fz3
Protein Names	Recommended name: Frizzled-3 Short name= dFz3
Expression Region	20-581
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.